

**MGM**  
I M P L A N T



**Implant**

*Plasma Activator*

## What is Implant *Plasma Activator*


**After the implant is activated after special surface activation, it is compatible with human tissue, which is conducive to cell attachment, and can accelerate bone tissue formation of hydrophilic implants to promote implant bone binding and shorten its healing cycle**





## USE COMPARISON CHART

### Hydrophilicity



The hydrophilicity was weak before plasma irradiation



Strong hydrophilicity after plasma irradiation



Strengthen adhesion and shorten healing time

### Powerfull Activation



## STORAGE

Implant surface is inevitably contaminated during cleaning, sterilization & storage, which decreases effective surface area

## PLASMA ACTIVATION

All contaminants are eliminated to regenerate effective surface area for delivery to patient under vacuum conditions



## COMPATIBILITY

All dental implant producters

